CE 0483

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# Instruction for use Monopolar Electrosurgical Electrode

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CE 0483

Copyright by Reger Medizintechnik GmbH  
All rights reserved  
Subject to technical changes without notice

Index: H

last update:  
2013-11-11

## ENGLISH

### 1 Scope

Monopolar electrosurgical electrodes Art. No.: 100-0 to 253-0, 255-0 to 691-0

Maximum rated voltage of accessory:

Article No.	U <sub>max</sub>
100-0 to 245-9 247-0 to 248-9 274-0 to 344-7 370-0 to 691-0	4,3 kVp
246-0 249-0 to 253-9 255-0 to 273-9 346-0 to 368-0	0,5 kVp

See also label or catalog.

In any combination with another electrosurgical accessory, the maximum rated voltage of the combination corresponds to the lowest rated voltage of the accessories used.

### 2 Authorized use

Only skilled medical personnel and personnel trained in this product are permitted to use the electrosurgical electrodes.

The electrosurgical electrodes are intended for open or endoscopic surgery and are used for cutting and coagulating biological tissue. They are not intended for direct contact with the heart or the central circulatory system. They are not intended for direct contact with the heart or central circulatory system. The electrosurgical electrodes are connected by an electrosurgical cable or an electrosurgical handle to the monopolar output of the electrosurgical generator. They are activated by means of a foot-operated switch or electrosurgical handle. It is recommended that a smoke evacuation system be used.

### 3 Safety notice

#### WARNING !

See these Instructions for Use, the label or the current product catalog for the maximum rated voltage of the electrosurgical electrodes.

If anything is unclear, contact the manufacturer.

Before initial use, new electrosurgical electrodes must be cleaned and also sterilized according to a validated procedure (DIN EN ISO 17665) (In this regard see Section 4 "Cleaning and sterilization"). A function test must be performed before each use (see Section 5 "Function test"). It must be ensured that the electrosurgical electrode is inserted firmly in the handle. This must be done carefully, in order to avoid damage to the electrosurgical electrode and/or injuries to the patient or surgical personnel. The electrode tip may be damaged if excessive force is applied. It is not permissible to activate the electrosurgical electrode as long as it is in contact with metal objects and/or optics. Throughout the complete procedure, care must be taken that no flammable substances are present in the immediate vicinity, since otherwise a danger of explosion exists. After the electrosurgical current has been turned off, the electrode tip may still be hot enough to cause burns.

The high-frequency current used in electrosurgery may interfere with cardiac pacemakers and implanted

heart defibrillators, and so affected patients must consult a cardiologist prior to the operation.

### 4 Cleaning and sterilization

In view of the design, the materials used and the intended use, a maximum limit cannot be defined for the number of cleaning and sterilization cycles that may be possible.

#### Preparation for cleaning:

Remove the electrosurgical electrodes from their packaging. Place them in a container/ device provided for cleaning/ sterilization.

It is not necessary to dismantle the electrosurgical electrodes.

#### Manual cleaning:

The electrosurgical electrodes must be disinfected immediately after each use. The electrosurgical electrodes should be thoroughly cleaned with a soft brush or synthetic fleece pad then rinsed, since otherwise particles or dried secretions may adhere to them. This could make subsequent cleaning and sterilization difficult or no longer possible. It must be ensured that areas with difficult access are thoroughly cleaned then rinsed several times. Highly alkaline cleaning agents with pH-levels above 10 and below 13 have no detrimental influence on the useful life.

#### Automated cleaning:

Only washing machines/ disinfectors with efficiency tested according to DIN EN ISO 15883 may be used. For cleaning, use only commercially available cleaners, which are suitable for cleaning of plastics parts. Thermal disinfection must be performed for at least 5 minutes at a minimum of 90°C [194°F] to a maximum to 95°C [203°F]. The manufacturer's instructions about the dosage of cleaning agent and selection of program must be followed. If chemical disinfection has been applied, a final rinse with distilled or demineralized water must be performed. The electrosurgical electrodes must be mounted securely and protected against mechanical damage during machine cleaning/ disinfection.

#### Sterilization:

The electrosurgical electrodes are designed exclusively for steam sterilization in autoclaves. The electrodes must be sterilized at a minimum of 134°C [273°F] and a maximum of 138°C [280°F] in saturated steam during a holding time of at least 5 minutes to at most 20 minutes, then dried in vacuum for at least 10 minutes. Sterilization must be performed in accordance with DIN EN ISO 17665 (Sterilization of Medical Devices in Moist Heat). It is not permissible to sterilize the electrodes with hot air, EO gas, gamma radiation or plasma.

**NOTE:** Before use, the electrosurgical electrodes must be cooled to room temperature.

#### Limitation of reconditioning:

The useful life of the product will depend on wear, damage and frequency of reconditioning.

### 5 Functional test

Before each use, the insulation of the electrosurgical-electrodes must be inspected for pressure points

or damage. It is not permissible to use electrosurgical electrodes exhibiting damage or pressure points.

### 6 Repair and modification

It is not permissible to repair defective electrosurgical electrodes. They must be replaced by new electrosurgical electrodes. Unauthorized modifications and repairs are strictly prohibited and will entail invalidation of the manufacturer's warranty.

### 7 Packaging, storage, transportation, handling

The electrosurgical electrodes must be stored in a clean and dry environment. They should be individually stored in a protective container with individual compartments or heat-sealed in film. The electrosurgical electrodes must always be handled with the utmost care during transportation, cleaning, upkeep, sterilization and storage. This applies in particular for fine tips and other sensitive areas. It is the operator's responsibility to ensure that the sterile condition is preserved after the sterilization process.

### 8 Returns

Returns will be accepted only if they are marked as "hygienically safe" or "not decontaminated" and have been securely packaged for shipping.

### 9 Disposal

The electrosurgical electrodes, the packaging material and the accessories must be disposed of in accordance with the regulations and laws specific to the country in which they are used.

### 10 About these Instruction for Use

Throughout the period of use of the electrosurgical electrodes, the Instructions for Use must be kept freely accessible for every user.

A current version may also be found at

[www.reger-med.de](http://www.reger-med.de).