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Monopolare HF-Elektrode (wiederverwendbar)

R001
Index: O

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Instructions for use

Monopolar Electrosurgical Electrode (reusable)

R001
Index: N

Last update:
2017-08-08

Copyright by Reger Medizintechnik GmbH
All rights reserved
Subject to technical changes without notice

ENGLISH

1 Scope

Monopolar electrosurgical electrodes Art. No.: 100-0 to 253-9, 255-0 to 691-0

Maximum rated voltage of accessory:

Article No.	U _{max}
100-0 to 245-9 247-0 to 248-10 272-0 to 344-7 350-2 373-0 to 691-0	4,3 kVp
246-0 249-0 to 253-9 255-0 to 268-0 345-0 to 350-1 351-0 to 368-0	0,5 kVp
270-0 to 271-10, 316-1 370-0, 371-0	1,3 kVp

See also label or catalogue.

In any combination with another electrosurgical accessory, the maximum rated voltage of the combination corresponds to the lowest rated voltage of the accessories used. (See also section 2 "Authorized Use").

2 Intended use

Only skilled medical professionals are permitted to use the electrosurgical electrodes. The electrosurgical electrodes are intended for open or endoscopic surgery and are used for cutting and coagulating biological tissue. They are not intended for direct contact with the heart or the central circulatory system. The electrosurgical electrodes are connected by a corresponding monopolar electrosurgical cable or a corresponding electrosurgical handle to the monopolar output of the electrosurgical generator. Monopolar electrosurgical electrodes are furnished either with a 1,6 mm, 2,4 mm 4 mm pin. 1,6 mm electrodes can be combined f.e. with REGER handle Art. No. 92098 or a comparable product. 2,4 mm electrodes can be combined f.e. with a REGER handle Art. No. 92251 or comparable product. 4 mm electrodes can be combined with f.e. REGER handle Art. No. 92247 or comparable product. The voltage of the used generator may not exceed the maximum rated voltage (See section 1 "Scope") of the electrosurgical electrodes and shall not exceed the frequency 4 MHz.

The electrodes are activated by means of a foot-operated switch or electrosurgical handle. It is recommended to use a smoke evacuation system. Monopolar electrosurgical electrodes can be operated with HF-generators of the following manufacturers:

▪ ERBE	▪ KLS Martin
▪ Covidien	▪ ValleyLab
▪ EMED	▪ Tekno
▪ BOWA	▪ Berchtold

Follow the instructions for use of the corresponding accessory and of the generator used.

3 Safety notice

WARNING !

See these Instructions for Use, the label or the current product catalog for the maximum rated voltage of the electrosurgical electrodes.

If anything is unclear, contact the manufacturer. Prior to each application, electrosurgical electrodes have to be cleaned, disinfected and sterilized according to a validated procedure (DIN EN ISO 17665) (In this regard see Section 4 "Cleaning, disinfection and sterilization"). A visual inspection must be performed before each use (see Section 5 "Visual inspection"). It must be ensured that the electrosurgical electrode is inserted firmly in the handle. This must be done carefully, in order to avoid damage to the electrosurgical electrode and/or injuries to the patient or surgical personnel. The electrode tip may be damaged if excessive force is applied. Prior to each usage, please verify if the designated neutral electrode is put on the patient correctly and whether the neutral electrode is connected with the HF-generator properly. It is not permissible to activate the electrosurgical electrode as long as it is in contact with metal objects and/or optics. Throughout the complete procedure, care must be taken that no flammable substances are present in the immediate vicinity, since otherwise a danger of explosion exists. After the electrosurgical current has been turned off, the electrode tip may still be hot enough to cause burns.

The high-frequency current used in electrosurgery may interfere with cardiac pacemakers and implanted heart defibrillators, and so affected patients must consult a cardiologist prior to the operation.

4 Cleaning, disinfection and sterilization

Reger Medizintechnik GmbH recommends the validated treatment processes described in the following. Equivalent deviant processes are possible. It is the sole responsibility of the user to safeguard the suitability of the actual applied procedure by suitable means, f.e. validation, routine examination, verification of material compatibility, etc. In view of the design, the materials used and the intended use of the product, a maximum limit cannot be defined for the number of cleaning, disinfection or sterilization cycles that may be possible. During proper use of the electrosurgical electrodes, the products are subject to natural wear and tear depending on type and duration of their application. Therefore a visual inspection has to be done prior to each application.

Preparation for cleaning:

Remove the electrosurgical electrodes from their packaging. Place them in a container/ device provided for cleaning/ sterilization. It is not necessary to dismantle the electrosurgical electrodes.

Manual pre-cleaning:

The electrosurgical electrodes must be disinfected immediately after each use. For pre-cleaning, please use water and where necessary aldehyde-free, non-fixing disinfectants. The electrosurgical electrodes should be thoroughly cleaned with a soft brush or synthetic fleece pad then rinsed, since otherwise particles or dried secretions may adhere to them. This could make subsequent cleaning and sterilization difficult or impossible. It must be ensured that areas with difficult access are thoroughly cleaned then rinsed several times.

Cavities and und lumen have to be rinsed intensively using at least 3 x 20 ml cold mains water (<40°C) with the aids of a rinsing adapter (f.e. from the company Medisafe), or with a syringe or with a water pressure pistol (>30 Sec.). This pre-cleaning step has to be done prior to the further manual cleaning or prior to cleaning with the cleaning and disinfection machine.

A) Manual cleaning and disinfection:

For cleaning, please prepare an immersion bath with a suitable fluid cleaning agent. Use a disinfection agent that is compatible to the cleaning agent, take care to use a disinfection agent that also can be used for immersion baths. Follow the instructions and recommendations of the manufacturer of the cleaning and disinfection agent. Only use agents that are suitable to clean medical devices made from metal and plastics with a pH-value of between 5.5 and 12.3. It is recommended to use the cleaning agent gigazyme® (Schülke & Mayr).

- According to the instructions of the manufacturer please prepare two separate immersion baths, one for cleaning and one for disinfection using the appropriate and respective agents.
- Put the HF electrode completely into the cleaner (f.e. 0.5% gigazyme®) in the ultrasonic bath applying an ultrasonic time of 5 Min. and a frequency of 35 kHz. Please follow the instructions of the manufacturer of the cleaning agent. Take care that the instruments do not touch any other parts inside the ultrasonic bath. Take care that sonic shadows are avoided inside the ultrasonic bath.
- Then, clean the HF electrodes with a soft brush under cold fluent town water (<40°C). Cavities and lumen have to be rinsed intensively by using a water pressure pistol (>30 Sec.) or similar. Afterwards rinse the HF electrodes thoroughly for at least 1 minute under mains water (<40°C) in order to remove any residues of the cleaning agent.
- Inspect the HF electrodes visually regarding remaining stain. In case soiling/stain is still detectable visually, repeat the aforementioned steps as long as it needs to remove any and all stain.
- Afterwards put the HF electrodes into the disinfection bath completely. Keep on the contact time set forth by the manufacturer of the disinfection agent, f.e. Korsolex Plus, 3% for 15 Min. Make sure that all areas and parts of the HF electrodes are covered by the disinfectant. Cavities and lumen have to be rinsed with the disinfectant for several times using a syringe f.e. and at least min. 3 times using 20 ml.
- Then, rinse the instruments thoroughly for at least 1 minute with demineralized water in order to remove any and all residues of the disinfectant.
- Dry the HF electrodes with a lint free wipe and sterile compressed air. Dry the cavities and lumen with sterile compressed air as well.

Highly alkaline cleaning agents with pH levels above 10 and below 13 have no detrimental influence on the lifetime of the electrosurgical electrodes. This pre-cleaning procedure has to be applied prior to further cleaning processes with cleaning and disinfecting machines.

Automated cleaning and disinfection:

Only apply cleaning and disinfecting machines that are efficiency tested according to DIN EN ISO 15883. Please refer to the information of the manufacturer of the cleaning or disinfecting agent, and only use agents that are suitable to be applied for medical devices made from metal and steel and have a pH-value between 5,5 and 12,3. It is recommended to use neodisher® mediclean forte (Dr. Weigert GmbH & Co. KG). Select the program for thermal disinfection. Regarding the program course, follow the instructions of use of the manufacturer of the agents. The electrosurgical electrodes must be mounted securely and protected against mechanical damage during machine cleaning/ disinfection Do not clean electrosurgical electrodes together with sharp-edged or pointed objects.

Put the electrosurgical electrodes in a suitable rinsing device.

Start the program course of the following features:

- Pre rinse with cold water for 1 Min.
- Emptying
- Pre rinse with cold water for 3 Min.
- Emptying
- Wash at 55°C for 5 Min. with a 0.5% alkaline cleaner
- Emptying
- Neutralisation with warm mains water(>40°C) and Neutralisator (0.1% Neodisher®Z) for 3 Min.
- Emptying
- Intermediate flusing: rinse with warm demineralised water (>40°C) for 2 Min.
- Emptying

Disinfection:

Proceed with the automated Thermic Disinfection according to the national requirements regarding the AO-value (see ISO 15883, AO-value >3000): 5 Min. for 92°C +/- 2°C.

Drying:

- 30 Min. at 90°C
- Remove rinsing adapter

At end of program course, remove the electrosurgical electrodes and verify if there are any staining residues. If staining/bonded residues are present, repeat the automated cleaning and disinfection process as long as no visible residues are present. Dry any cavities and insufficiently dried areas with sterile compressed air < 2 bar. After the electrosurgical electrodes are taken out (after additional re-drying on a clean place), immediately pack them in a single-use sterilization package (double package) made of paper/foil or put them in a sterilization container (according to DIN EN ISO 11607 and DIN EN 868).

Sterilization:

It is only allowed to sterilize instruments that have been cleaned and disinfected beforehand. Electrosurgical electrodes are designed exclusively for steam sterilization in autoclaves (fractionated fore-vacuum procedure with sufficient product drying). Electrosurgical electrodes have to be sterilized at a minimum of 134°C and maximum of 137°C in saturated steam during a holding time of at least 5 minutes to at most 20 minutes, then dried in a vacuum and sterilized for at least 10 minutes. Sterilizer according to valid national Norms and Guidelines (f.e. DIN EN 13060 or DIN EN 285). The sterilization process has been validated according to DIN EN ISO 17665. Please follow the recommendations and instructions of the sterilizer manufacturer regarding loading, handling and drying times. It is not permissible to sterilize the electrosurgical electrodes with hot air, EO gas, gamma radiation or plasma.

NOTE: Before use, the electrosurgical electrodes must be cooled to room temperature. It is the sole responsibility of the user to safeguard and maintain the sterile condition of the electrosurgical electrodes after the sterilization process.

In case the aforementioned described and recommended chemicals and machines for manual pre-cleaning and/or automated cleaning and disinfection are not available, it is the responsibility of the user to validate their procedure. Also, if another sterilization procedure is chosen, the procedure deviant from the procedure described in paragraph 4 has to be validated by the user accordingly.

Limitation of reconditioning:

The lifetime of the product will depend on wear, damage and frequency of reconditioning.

5 Visual inspection

Before each use, the insulation of the electrosurgical electrodes must be inspected for pressure points or damage. It is not permissible to use electrosurgical electrodes exhibiting damage or pressure points.

6 Repair and modification

It is not permissible to repair defective electrosurgical electrodes. They must be replaced by new electrosurgical electrodes. Unauthorized modifications and repairs are strictly prohibited and will entail invalidation of the manufacturer's warranty.

7 Packaging, storage, transportation, handling

The electrosurgical electrodes must be stored in a clean and dry environment. They should be individually stored in a protective container with individual compartments or heat-sealed in film. The electrosurgical electrodes must always be handled with the utmost care during transportation, cleaning, disinfection, upkeep, sterilization and storage. This applies in particular for fine tips and other sensitive areas. It is the operator's responsibility to ensure that the sterile condition is preserved after the sterilization process.

8 Returns

Returns will be accepted only if they are marked as "hygienically safe" or "not contaminated" and have been securely packaged for shipping.

9 Disposal

The electrosurgical electrodes, the packaging material and the accessories must be disposed of in accordance with the regulations and laws specific to the country in which they are used.

10 About this Instructions for Use

Throughout the period of use of the electrosurgical electrodes, the Instructions for Use must be kept freely accessible for every user. A current version may also be found on www.reger-med.de.

©2016 Reger Medizintechnik GmbH. All rights reserved. All company names, trademarks, service marks, trade names, product names and logos appearing on this document are the property of their respective owners and should be treated accordingly.