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# Instructions for Use Electrosurgical Cables and Adapters

## ENGLISH

### 1 Scope

**Bipolar HF-cables/adapters Art.-No.:** 92015 to 92020, 92023 to 92026, 92029 to 92032, 92033, 92037, 92038, 92065 to 92068, 92073 to 92075, 92084 to 92088, 92092, 92094, 92110, 92111, 92113, 92115, 92116, 92127, 92128, 92130 to 92133, 92136 to 92138, 92140 to 92142, 92152, 92154, 92155, 92185, 92187, 92190, 92192, 92196 to 92198, 92204 to 92212, 92224, 92225, 92230, 92232, 92389 to 92391, 92394, 92397, 92398, 92451 to 92453, 92498, 92504

**Monopolar HF-cables/adapters Art.-No.:** 92002 to 92014, 92039 to 92041, 92050 to 92052, 92061 to 92063, 92070 to 92072, 92079 to 92083, 92089, 92090, 92092, 92096 to 92098, 92100, 92105 to 92107, 92109, 92112, 92119, 92161 to 92164, 92175 to 92178, 92180 to 92182, 92184, 92191, 92195, 92199, 92201, 92203, 92220, 92222, 92223, 92226, 92395, 92396, 92399, 92901

**Monopolar HF-cables for neutral electrode, Art.-No.:** 92076 to 92078, 92108, 92109, 92139, 92145, 92179, 92344, 92345, 92497, 92499.

Maximum rated voltage of accessory:

Designation	U <sub>max</sub>
Monopolar cable/adaptor	4,3 kVp
Bipolar cable/adaptor	0,8 kVp
Neutral electrode cable	4,3 kVp

See also label or catalog.

In any combination with another electrosurgical accessory, the maximum rated voltage of the combination corresponds to the lowest rated voltage of the accessories used.

### 2 Intended use

Only skilled medical personnel are permitted to use the electrosurgical cables/ adapters.

Electrosurgical cables/ adapters are used to connect electrosurgical-accessories for monopolar or bipolar application to an electrosurgical generator.

Follow the instructions for use of the corresponding electrosurgical accessory and of the generators used.

### 3 Safety notice

#### WARNING !

See these Instructions for Use, the label or the current product catalog for the maximum rated voltage of the accessory.

If anything is unclear, contact the manufacturer. Before initial use, new electrosurgical cables/ adapters must be cleaned, disinfected and sterilized according to a validated procedure (DIN EN 17665) (In this regard see Section 4 "Cleaning, disinfection and sterilization"). A function test and a visual inspection must be performed before each use (see Section 5 "Visual inspection and function test"). If must be ensured that the electrosurgical cable/ adapter is connected firmly and correctly to the generator. This needs to be done carefully in order to avoid damages on HF-cables or adapters. Throughout the complete procedure, care must be taken that no flammable substances are present in the immediate vicinity, since otherwise a danger of explosion exists.

The high-frequency-current used in electrosurgery may interfere with cardiac pacemakers and implanted heart defibrillators, and so affected patients must consult a cardiologist prior to the operation.

### 4 Cleaning, disinfection and sterilization

Reger Medizintechnik GmbH recommends the validated treatment processes described in the following. Equivalent deviant processes are possible. It is the sole responsibility of the user to safeguard the suitability of the actual applied procedure by suitable means, i.e. validation, routine examination, verification of material compatibility, etc.

In view of the design, the materials used and the intended use of the product, a maximum limit cannot be defined for the number of cleaning, disinfection or sterilization cycles that may be possible. During proper use of the electrosurgical cables/adapters, the products are subject to natural wear and tear depending on type

and duration of their application. Therefore a sight check and functional test has to be done prior to each application.

#### Preparation for cleaning:

Where applicable, remove adapter from the cable. Unpack cables/adapters and place them in a compartment suitable for pre-cleaning.

#### Manual pre-cleaning:

The electrosurgical cables/adapters must be disinfected immediately after each use. For pre-cleaning, use water and if necessary, aldehyde-free, non-fixing disinfectants. The cables/adapters should be thoroughly cleaned with a soft brush or synthetic fleece pad then rinsed, since otherwise particles or dried secretions may adhere to them. This could make subsequent cleaning and sterilization difficult or impossible. It must be ensured that areas with difficult access are thoroughly cleaned then rinsed several times. Highly alkaline cleaning agents with pH levels above 10 and below 13 have no detrimental influence on the lifetime of the electrosurgical cables/adapters. This pre-cleaning procedure has to be applied prior to further cleaning processes with cleaning and disinfecting machines.

#### Automated cleaning and disinfection:

Only apply cleaning and disinfecting machines that are efficiency tested according to DIN EN ISO 15883. Please refer to the information of the manufacturer of the cleaning or disinfecting agent, and only use agents that are suitable to be applied for medical devices made from metal and steel and have a pH-value between 5,5 and 12,3. It is recommended to use neodisher® mediclean forte (Dr. Weigert GmbH & Co. KG). Do not use neutralization agents. Select the program for thermal disinfection. Regarding the program course, follow the instructions of use of the manufacturer of the agents. Do not clean electrosurgical handles together with sharp-edged or pointed objects.

Put the cables/adapters in a suitable rinsing device and take care that the cable is not folded or squeezed during cleaning and disinfection. Start the program course of the following features:

- Thermal disinfection: 5 to 10 minutes at 90°C to 93 °C (tolerance according to DIN EN ISO 15883-1 or -2), AO ≥ 3000.
- Final rinsing with demineralized water.
- Drying of product

At end of program course, remove the handles and verify if there are any staining residues. If staining/bonded residues are present, repeat the automated cleaning and disinfecting process as long as no visible residues are present. Dry any cavities and insufficiently dried areas with sterile compressed air < 2 bar. After the cables/adapters are taken out (after additional re-drying on a clean place), immediately pack them in a single-use sterilization package (single or double package) made of paper/foil or put them in a sterilization container (according to DIN EN ISO 11607 and DIN EN 868).

#### Sterilization:

It is only allowed to sterilize cables/adapters that have been cleaned and disinfected beforehand. Electrosurgical cables/adapters are designed exclusively for steam sterilization in autoclaves (fractionated fore-vacuum procedure with sufficient product drying). HF-cables/adapters have to be sterilized at a minimum of 132°C and maximum of 138°C in saturated steam during a holding time of at least 3 minutes to at most 20 minutes, then dried in a vacuum and sterilized for at least 10 minutes. Select a sterilization device according to valid national Norms and Guidelines (i.e. DIN EN 13060 or DIN EN 285). The sterilization process has been validated according to DIN EN ISO 17665. Please follow the recommendations and instructions of the sterilization device manufacturer regarding loading, handling and drying times. It is not permissible to sterilize the electrosurgical cables/adapters with hot air, EO gas, gamma radiation or plasma.

**NOTE:** Before use, the electrosurgical cables/adapters must be cooled to room temperature. It is the sole responsibility of the user to safeguard and maintain the sterile condition of the cables/adapters after the sterilization process. In case the aforementioned described and recommended chemicals and machines for manual pre-cleaning and/or automated cleaning and disinfection are not available, it is the responsibility of the user to validate their

procedure. Also, if another sterilization procedure is chosen, the procedure deviant from the procedure described in paragraph 4 has to be validated by the user accordingly.

**Limitation of reconditioning:** The lifetime of the product will depend on natural wear and tear, damage and frequency of reconditioning.

### 5 Visual and function test

Prior to each use, the insulation of the electrosurgical cables/ adapters must be inspected for pressure points or damage. In addition, an electrical continuity test must be performed before each use. It is not permissible to use damaged electrosurgical cables/ adapters or electrosurgical cables/ adapters exhibiting pressure points.

### 6 Repair and modification

It is not permissible to repair defective electrosurgical cables/ adapters. They must be replaced by new electrosurgical-cables/ adapters. Unauthorized modifications and repairs are strictly prohibited and will entail invalidation of the manufacturer's warranty.

### 7 Packaging, storage, transportation, handling

The electrosurgical cables/ adapters must be stored in a clean and dry environment. They should be individually stored in a protective container with individual compartments or heat-sealed in film. The electrosurgical-cables/ adapters must always be handled with the utmost care during transportation, cleaning, upkeep, sterilization and storage. This applies in particular for fine tips and other sensitive areas. It is the operator's responsibility to ensure that the sterile condition is preserved after the sterilization process.

### 8 Returns

Returns will be accepted only if they are marked as "hygienically safe" or "not contaminated" and have been securely packaged for shipping.

### 9 Disposal

The electrosurgical-cables/ adapters, the packaging material and the accessories must be disposed of in accordance with the regulations and laws specific to the country in which they are used.

### 10 About this Instructions for Use

Throughout the period of use of the Instructions for Use of the electrosurgical cables/adapters must be kept freely accessible for every user.

A current version may also be found on

[www.reger-med.de](http://www.reger-med.de).

©2018.REGER Medizintechnik GmbH. All rights reserved. All company names, trademarks, service marks, trade names, product names and logos appearing on this document are the property of their respective owners and should be treated accordingly.