DEUTSCH

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Monopolare Laparoskopie-Elektroden Art.-Nr.:

700-0 bis 769-0 und 784-0

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
Copyright by Reger Medizintechnik GmbH
All rights reserved
Subject to technical changes without notice

ENGLISH

1 Scope

Monopolar laparoscopic electrodes Art. No.:

700-0 to 769-0 and 784-0

Maximum rated voltage of accessory:

Article No.	U _{max}
700-0 to 769-0 784-0	4.3 kVp

See also label or catalog.

In any combination with another electrosurgical accessory, the maximum rated voltage of the combination corresponds to the lowest rated voltage of the accessories used. (See also section 2 "Authorized Use")

2 Intended use

Only skilled medical personnel are permitted to use the laparoscopic electrodes. The laparoscopic electrodes are intended for laparoscopic use. They are designed for cutting and coagulating biological tissue. They are not intended for direct contact with the heart or central circulatory system. The laparoscopic electrodes are connected by an electrosurgical cable to the monopolar output of the electrosurgical generator. Laparoscopic electrodes are furnished with a 4 mm pin and can be combined with REGER handle Art. No. 92247 f.e. or with a comparable product. Laparoscopic electrodes may be operated only with HF-generators that are furnished with a monopolar exit and with a max. frequency of 4 MHz. Exceptions are the laparoscopic electrodes Art.-No. 722-0, 725-0 and 726-0. These laparoscopic electrodes are not suitable for RF-generators, whereas the frequency of 400 kHz shall not be exceeded. Laparoscopic electrodes are activated by a foot-operated switch. Laparoscopic electrodes can be operated with HF-generators of the following manufacturers:

ERBE	KLS Martin
Covidien	ValleyLab
EMED	Tekno
BOWA	Berchtold

Follow the instructions for use of the corresponding accessory and of the generator used.

3 Safety notice !WARNING!

See these Instructions for Use, the label or the current product catalog for the maximum rated voltage of the laparoscopic electrodes.

If anything is unclear, contact the manufacturer.

Prior to each application, laparoscopic electrodes have to be cleaned, disinfected and sterilized (DIN EN ISO 17665) according to a validated procedure, (refer also to paragraph 4 „Cleaning disinfection and sterilization“). Prior to each application, a visual inspection has to be done (see paragraph 5 „Visual Inspection“). It must be ensured that the laparoscopic electrode is firmly connected to the electrosurgical cable. This must be done carefully, in order to avoid damage to the laparoscopic electrode and/ or injuries to the patient or surgical personnel. The electrode tip may be damaged if excessive force is applied. Prior to each usage, please verify if the designated neutral electrode is put on the patient correctly and whether the neutral electrode is connected with the HF-generator properly. It is not permissible to activate the laparoscopic electrode as long as it is in contact with metal objects and/ or optics. Throughout the complete procedure, care must be taken that no flammable substances are present in the immediate vicinity, since otherwise a danger of explosion exists. After the electrosurgical current has been turned off, the electrode tip may still be hot enough to cause burns. The high-frequency current used in electrosurgery may interfere with cardiac pacemakers and implanted heart defibrillators, and so affected patients must consult a cardiologist prior to the operation.

4 Cleaning, disinfection and sterilization

Reger Medizintechnik GmbH recommends the validated treatment processes described in the following. Equivalent deviant processes are possible. It is the sole responsibility of the user to safeguard the suitability of the actual applied procedure by suitable means, f.e. validation, routine examination, verification of material compatibility, etc.

In view of the design, the materials used and the intended use of the product, a maximum limit cannot be defined for the number of cleaning, disinfection or sterilization cycles that may be possible. During proper use of the laparoscopic electrodes, the products are subject to natural wear and tear depending on type and duration of their application. Therefore a visual inspection has to be done prior to each application.

Preparation for cleaning:

Remove the laparoscopic electrodes from their packaging. Place them in a container/ device provided for cleaning/ sterilization.

It is not necessary to dismantle the laparoscopic electrodes.

Manual pre-cleaning:

The laparoscopic electrodes must be disinfected immediately after each use. For pre cleaning use water and aldehyde-free, non-fixing disinfectant if necessary. The laparoscopic electrodes should be thoroughly cleaned with a soft brush or synthetic fleece pad then rinsed, since otherwise particles or dried secretions may adhere to them. This could make subsequent cleaning and sterilization difficult or impossible. It must be ensured that areas with difficult access are thoroughly cleaned then rinsed several times.

Cavities and und lumen (f. e. rinsing/suction channel) have to be rinsed intensively using at least 3 x 20 ml cold mains water (<40°C) with the aids of

a rinsing adapter (f.e. from the company Medisafe), or with a syringe or with a water pressure pistol (>30 Sec.).

This pre-cleaning step has to be done prior to the further manual cleaning or prior to cleaning with the cleaning and disinfection machine.

A) Manual cleaning and disinfection:

For cleaning please prepare an immersion bath with a suitable fluid cleaning agent. Use a disinfection agent that is compatible to the cleaning agent. Take care to use a disinfection agent that also can be used for immersion baths. Follow the instructions and recommendations of the manufacturer of the cleaning and disinfection agent. Only use agents that are suitable to clean medical devices made from metal and plastics with a pH-value of between 5.5 and 12.3. It is recommended to use the cleaning agent gigazyme® (Schülke & Mayr).

1. According to the instructions of the manufacturer please prepare two separate immersion baths, one for cleaning and one for disinfection using the appropriate and respective agents.
2. Put the laparoscopic electrode completely into the cleaner (f.e. 0.5% gigazyme®) in the ultrasonic bath applying an ultrasonic time of 5 Min. and a frequency of 35 kHz. Please follow the instructions of the manufacturer of the cleaning agent. Take care that the instruments do not touch any other parts inside the ultrasonic bath. Take care that sonic shadows are avoided inside the ultrasonic bath.
3. Then, clean the laparoscopic electrodes with a soft brush under cold fluent town water (<40°C). Cavities and lumen have to be rinsed intensively by using a water pressure pistol (>30 Sec.) or similar. Afterwards rinse the laparoscopic electrodes thoroughly for at least 1 minute under mains water (<40°C) in order to remove any residues of the cleaning agent.
4. Inspect the laparoscopic electrodes visually regarding remaining stain. In case soiling/stain is still detectable visually, repeat the aforementioned steps as long as it needs to remove any and all stain.
5. Afterwards put the laparoscopic electrodes into the disinfection bath completely. Keep on the contact time set forth by the manufacturer of the disinfection agent, f.e. Korsolox Plus, 3% for 15 Min. Make sure that all areas and parts of the laparoscopic electrodes are covered by the disinfectant. Cavities and lumen have to be rinsed with the disinfectant for several times using a syringe f.e. and at least min. 3 times using 20 ml.
6. Then, rinse the instruments thoroughly for at least 1 minute with demineralized water in order to remove any and all residues of the disinfectant. Regarding laparoscopic electrodes that have a rinsing/suction channel, please rinse this channel at least 3 times using each 20 ml demineralized water. Rinse the channel with a rinsing adapter (f.e. from the company Medisafe) or with a syringe.
7. Dry the laparoscopic electrodes with a lint free wipe and sterile compressed air. Dry the cavities and channels with sterile compressed air as well.

B) Automated cleaning and disinfection:

Only apply cleaning and disinfecting machines that are efficiency tested according to DIN EN ISO 15883. Please refer to the information of the manufacturer of the cleaning or disinfecting agent, and only use agents that are suitable to be applied for medical devices made from metal and steel and have a pH-value between 5.5 and 12.3. It is recommended to use neodisher® mediclean forte (Dr. Weigert GmbH & Co. KG). Select the program for thermal disinfection. Regarding the program course, follow the instructions of use of the manufacturer of the agents. Do not clean laparoscopic electrodes together with sharp-edged or pointed objects.

Put the laparoscopic electrodes in a suitable rinsing device. Connect the instruments to the rinsing device using a rinsing adapter (f.e. from the company Medisafe).

- Pre rinse with cold water for 1 Min.
- Emptying
- Pre rinse with cold water for 3 Min.
- Emptying
- Wash at 55°C for 5 Min. with a 0.5% alkaline cleaner
- Emptying
- Neutralisation with warm mains water(>40°C) and Neutralisator (0.1% Neodisher®Z) for 3 Min.
- Emptying
- Intermediate flusing: rinse with warm demineralised water (>40°C) for 2 Min.
- Emptying

Disinfection:

Proceed with the automated Themic Disinfection according to the national requirements regarding the AO-value (see ISO 15883, AO-value >3000): 5 Min. for 92°C +/- 2°C.

Drying:

- 30 Min. at 90°C
- Remove rinsing adapter

At the end of the program course, remove the laparoscopic electrodes and verify if there are any staining residues. If staining/bonded residues are present, repeat the automated cleaning and disinfecting process as long as no visible residues are present. Dry any cavities and insufficiently dried areas with sterile compressed air < 2 bar. After the laparoscopic electrodes are taken out (after additional re-drying on a clean place), immediately pack

them in a single-use sterilization package (double package) made of paper/foil or put them in a sterilization container (according to DIN EN ISO 11607 and DIN EN 868).

Sterilization:

It is only allowed to sterilize instruments that have been cleaned and disinfected beforehand. Laparoscopic electrodes are designed exclusively for steam sterilization in autoclaves (fractionated fore-vacuum procedure with sufficient product drying). Laparoscopic electrodes have to be sterilized at a minimum of 134°C and maximum of 137°C in saturated steam during a holding time of at least 3 minutes to at most 20 minutes, then dried in a vacuum and sterilized for at least 10 minutes. Sterilization according to valid national Norms and Guidelines (f.e. DIN EN 13060 or DIN EN 285). The sterilization process has been validated according to DIN EN ISO 17665. Please follow the recommendations and instructions of the sterilization manufacturer regarding loading, handling and drying times. It is not permissible to sterilize the laparoscopic electrodes with hot air, EO gas, gamma radiation or plasma.

NOTE: Before use, laparoscopic electrodes must be cooled to room temperature. It is the sole responsibility of the user to safeguard and maintain the sterile condition of the laparoscopic electrodes after the sterilization process. In case the aforementioned described and recommended chemicals and machines for manual pre-cleaning and/or automated cleaning and disinfection are not available, it is the responsibility of the user to validate their procedure. Also, if another sterilization procedure is chosen, the procedure deviant from the procedure described in paragraph 4 has to be validated by the user accordingly.

Limitation of reconditioning:

In view of the design, the materials used and the intended use of the product, a maximum limit cannot be defined for the number of cleaning, disinfection or sterilization cycles that may be possible. The lifetime of the product will depend on natural wear and tear, damage and frequency of reconditioning.

5 Visual inspection

Before each use, the insulation of the laparoscopic electrodes must be inspected for pressure points or damage. It is not permissible to use laparoscopic electrodes exhibiting damage or pressure points.

6 Repair and modification

It is not permissible to repair defective laparoscopic electrodes. They must be replaced by new laparoscopic electrodes. Unauthorized modifications and repairs are strictly prohibited and will entail invalidation of the manufacturer's warranty.

7 Packaging, storage, transportation, handling

The laparoscopic electrodes must be stored in a clean and dry environment. They should be individually stored in a protective container with individual compartments or heat-sealed in film. The laparoscopic electrodes must always be handled with the utmost care during transportation, cleaning, upkeep, sterilization and storage. This applies in particular for fine tips and other sensitive areas. It is the operator's responsibility to ensure that the sterile condition is preserved after the sterilization process.

8 Returns

Returns will be accepted only if they are marked as "hygienically safe" or "not contaminated" and have been securely packaged for shipping.

9 Disposal

The laparoscopic electrodes, the packaging material and the accessories must be disposed of in accordance with the regulations and laws specific to the country in which they are used.

10 About this Instructions for Use

Throughout the period of use of the laparoscopic electrodes, the Instructions for Use must be kept freely accessible for every user.

A current version may also be found at

www.reger-med.de.

©2017 Reger Medizintechnik GmbH. All rights reserved. All company names, trademarks, service marks, trade names, product names and logos appearing on this document are the property of their respective owners and should be treated accordingly.