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Uterus-Manipulator Art. Nr. 75001

R014
Index: F

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Copyright by Reger Medizintechnik GmbH
All rights reserved
Subject to technical modifications without further notice

ENGLISH

1 Scope

Uterus Manipulator Art.-No.: 75001

Surgical instrument for application of a liquid. When combining the instrument with any accessory the safety notices of the several accessory components have to be observed as well.

2 Intended Use

Only skilled medical professionals are allowed to use the Uterus Manipulator.

The Uterus Manipulator serves for the positioning of the uterus during laparoscopic applications in the underbelly and for the introduction of a liquid in the uterus and tubes in order to test their patency. Only commercially available medical single use syringes that are furnished with a Luer-Lock connector are allowed to be applied.

3 Safety Notices

WARNING !

The intended use of the instrument can be retrieved from this Instructions for Use. If anything is unclear contact the manufacturer.

The Uterus Manipulator is a single use instrument. Prior to usage a new Uterus Manipulator has to be cleaned, disinfected and sterilized according to a validated procedure (DIN EN ISO 17665) (see also paragraph 4 "Cleaning, disinfection and sterilization"). Prior to usage a visual inspection and a function test has to be done (see paragraph 5 "visual inspection and function test"). It has to be safeguarded that all connections are combined firmly and correctly with each other. This has to be done carefully in order to avoid damages on the Uterus Manipulator, on the accessory and to avoid injuries of the patient and/or medical personnel.

Do not apply the Uterus Manipulator in case the patient suffers from heavy bleeding, active infection of the pelvis and/or existing pregnancy.

4 Cleaning, disinfection and sterilization

The Uterus Manipulator only is allowed to be used one time. It is a Single Use device. Therefore a reesterilization after usage is strictly prohibited.

REGER Medizintechnik GmbH recommends an automated cleaning and disinfection. Due to minor efficacy a manual cleaning is not recommended. Do not use sharp edged objects to clean the instrument. After disinfection the disinfectant agent has to be removed thoroughly by rinsing.

Reger Medizintechnik GmbH recommends the validated treatment processes described in the following. Equivalent deviant processes are possible. It is the sole responsibility of the user to safeguard the suitability of the actual applied procedure by suitable means, f.e. validation, routine examination, verification of material compatibility, etc.

Preparation for cleaning:

Prior to its usage the Uterus Manipulator has to be taken out of its package. It has to be placed in a container for cleaning/sterilization. During cleaning/disinfection the instrument has to be stored safely in order to protect it from mechanical damages.

It is not necessary to dismantle the Uterus Manipulator.

Manual pre-cleaning :

The exterior of the instrument has to be cleaned with a moist lint free wipe. During this cleaning the plastic cone has to be moved back and forth in axial direction several times.

The interior channel has to be rinsed with at least x 20 ml mains water using a rinsing adapter (f.e. from the company Medisafe) via the Luer-Lock connector and using a standard syringe without cannula. This pre-cleaning step has to be done prior to the further automated cleaning with the cleaning and disinfection machine.

Automated cleaning and disinfection:

Only apply cleaning and disinfecting machines that are efficiency tested according to DIN EN ISO 15883. Please refer to the information of the manufacturer of the cleaning or disinfecting agent, and only use agents that are suitable to be applied for medical devices made from metal and steel and have a pH-value between 5,5 and 12,3. It is recommended to use neodisher® mediclean forte (Dr. Weigert GmbH & Co. KG). Select the program for thermal disinfection. Regarding the program course, follow the instructions of use of the manufacturer of the agents. The electrosurgical electrodes must be mounted securely and protected against mechanical damage during machine cleaning/ disinfection Do not clean electrosurgical electrodes together with sharp-edged or pointed objects.

Put the instrument in a suitable rinsing device. Use a rinsing adapter (f.e. from the company Medisafe) to connect the instrument with the rinsing device.

Start the program course with following features:

- 1 min. pre-rinsing with cold water
- Emptying
- 3 min. pre-rinsing with cold water
- Emptying
- 5 min washing at 55° C with 0.5 % alkaline cleaner
- Emptying
- 3 min neutralisation with warm mains water (>40°C) and neutralisator (0.1% Neodisher®Z)
- Emptying
- 2 min rinsing with warm demineralized water (>40°C)
- Emptying

Disinfection :

Proceed with the automated Thermic Disinfection according to the national requirements regarding the AO-value (see ISO 15883, AO-value >3000): 5 Min. for 92°C +/- 2°C.

Drying:

- 30 Min. at 90°C
- Remove rinsing adapter

At end of program course, remove the electrosurgical electrodes and verify if there are any staining residues. If staining/bonded residues are present, repeat the automated cleaning and disinfecting process as long as no visible residues are present. Dry any cavities and insufficiently dried areas with sterile compressed air < 2 bar. After the electrosurgical electrodes are taken out (after additional re-drying on a clean place), immediately pack them in a single-use sterilization package (double package) made of paper/foil or put them in a sterilization container (according to DIN EN ISO 11607 and DIN EN 868).

Sterilization:

Proceed with steam sterilization in the autoclave for at minimum 134° C and maximum 138° C in saturated steam. Holding time min. 5 min. to at least 20 min. Drying in the vacuum for min. 10 Min. DIN EN ISO 17665 (sterilization of medical devices in humid heat) has to be observed. Do not sterilize in hot air. The Uterus Manipulator is not allowed to be sterilized with hot air, EO gas, gamma radiation or plasma.

Remark: Prior to usage the Uterus Manipulator has to cool down to room temperature.

Limitation of reconditioning:

Due to the chosen design, the materials used and the intended use of the instrument the Uterus Manipulator only is allowed to be applied one time.

5 Visual inspection and function test

Prior to usage the instrument has to be checked for damages. Furthermore a function test has to be done. During the function test it has to be checked whether the cohesion of the black plastic cone is sufficient. To do so, place the cone over the length marks. Screw on completely the union nut on the thread of the cone. Try to move the cone in axial direction. It should only be possible that the cone slips to the next length mark at the maximum, not further. It is not allowed to use damage instruments and/or instruments on which the cone is not closing correctly.

6 Repair and Modification

It is not permissible to repair defective Uterus Manipulators. They must be replaced by new Uterus Manipulators. Unauthorized modifications and repairs are strictly prohibited and will entail invalidation of the manufacturer's warranty.

7 Packing, Storage, Transport and Handling

The electrosurgical electrodes must be stored in a clean and dry environment. They should be individually stored in a protective container with individual compartments or heat-sealed in film. The electrosurgical electrodes must always be handled with the utmost care during transportation, cleaning, disinfection, upkeep, sterilization and storage. This applies in particular for fine tips and other sensitive areas. It is the operator's responsibility to ensure that the sterile condition is preserved after the sterilization process.

8 Returns

Returns will be accepted only if they are marked as "hygienically safe" or "not contaminated" and have been securely packaged for shipping.

9 Disposal

The electrosurgical electrodes, the packaging material and the accessories must be disposed of in accordance with the regulations and laws specific to the country in which they are used.

10 About this Instructions for Use

Throughout the period of use of the electrosurgical electrodes, the Instructions for Use must be kept freely accessible for every user.

A current version may also be found at

www.reger-med

©2016.Reger Medizintechnik GmbH. All rights reserved. All company names, trademarks, service marks, trade names, product names and logos appearing on this document are the property of their respective owners and should be treated accordingly.